



Cylinder, graduated, 1000 mL, w/ handles

Ürün Numarası: 2612953
TRY Fiyat: Bize Ulaşın
Mevcut

One-piece polypropylene construction with a molded in handle top ensure a secure, solid grip and superior control for safer, more convenient pouring. The “Total Contained, Total Delivered ”double scale is accurate and easy to read (cylinder will not create a meniscus), while the hexagonal base prevents tipping or rolling. Hand opening is 89x 29 mm (3 x 11/8"). Steam autoclavable upto 121°C (250°F).

Teknik özellikler

Gövde malzemesi:	Polypropylene
Hacim / Paket boyutu:	1 each
Kapasite:	1000 mL
Ölçüm metodu:	To Deliver/To Contain
Sertifikalar:	None
Tür:	Çift Metrik