



## Cylinder, graduated, 25 mL +- 0.3 mL, 0.5 mL divisions (white markings)

Ürün Numarası: 50840  
TRY Fiyat: Bize Ulaşın  
Mevcut

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes  $\geq 25$  mL.

### Teknik özellikler

|                       |            |
|-----------------------|------------|
| Gövde malzemesi:      | Glass      |
| Hacim / Paket boyutu: | 1 each     |
| Kapasite:             | 25 mL      |
| Ölçüm metodu:         | To Contain |
| Sertifikalar:         | Class B    |
| Tür:                  | Tek Metrik |