



Cylinder, graduated, 10 mL +- 0.1 mL, 0.2 mL divisions (white markings)

Ürün Numarası: 50838
TRY Fiyat: Bize Ulaşın
Mevcut

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes ≥ 25 mL.

Teknik özellikler

Gövde malzemesi:	Glass
Hacim / Paket boyutu:	1 each
Kapasite:	10 mL
Ölçüm metodu:	To Deliver
Sertifikalar:	Class B
Tür:	Tek Metrik
Yükseklik:	135 mm