



## Bottle, BOD, 300 mL, 24/pk, serialized: 1-24

Ürün Numarası: 2898700  
TRY Fiyat: Bize Ulaşın  
Mevcut

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included Serialized bottles have silk-screened, sequential numbers (1-24) and include a bar code.

### Teknik özellikler

Body Material:	Glass
Hacim / Paket boyutu:	24 /pk
Kapasite:	300 mL
Malzeme:	Borosilicate Glass
Package size UOM:	24/pk, Serialized: 1 - 24