



Cylinder, graduated, polymethylpentene, 25 mL ± 0.3 mL, 0.5 mL divisions

Ürün Numarası: 217240
TRY Fiyat: Bize Ulaşın
2 hafta içinde teslimat

This cylinder features transparent, nontoxic, noncontaminating, and biologically inert material with excellent chemical resistance. Meets ASTM E1272 Class B accuracy requirements, ISO Standard 6706, and suitable for food and beverage use (PMP resin meets requirements of CFR 21, Part 177.1520). Calibrated “To Contain/To Deliver” at 20°C (68°F). Shatterproof and thermally resistant to 135°C (275°F). Absence of meniscus for easy-to-read levels on front-to-back raised graduations.

Teknik özellikler

Gövde malzemesi:	Polymethylpentene
Hacim / Paket boyutu:	1 each
Kapasite:	25 mL
Ölçüm metodu:	To Deliver/To Contain
Sertifikalar:	Class B
Tür:	Tek Metrik