



Bottle, BOD, 60 mL capacity, marked 30 mL graduation line

Ürün Numarası:	190902
TRY Fiyat:	Bize Ulaşın
Mevcut	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 mL sizes).

Teknik özellikler

Body Material:	Borosilicate Glass
Hacim / Paket boyutu:	1 each
Kapasite:	60 mL
Malzeme:	Borosilicate Glass
Package size UOM:	1pk