



Bottle, wash, diamond alert, 500 mL, Variety 6/pk

Ürün Numarası:

2613049

TRY Fiyat:
Mevcut

Bize Ulaşın

Wide-mouthed and color-coded wash bottles provide safe and positive identification of common laboratory solutions. Standard Diamond Alert symbols indicate type of hazard contents may represent. Polyethylene bottles with color-coded polypropylene screw caps and spray-tip assembly. Variety 6-pack includes one bottle each for water, acetone, methanol, isopropanol, toluene, and ethanol.

Teknik özellikler

| | |
|-----------------------|---------------|
| Birim: | 6 pk |
| Hacim / Paket boyutu: | 6 pk |
| Kapasite: | 500 mL |
| Malzeme: | LDPE |
| Renk: | Transparent |
| Stem: | Blue Cap/Stem |