



Cylinder, graduated, mixing, certified, 25 mL

Ürün Numarası: 2636340
TRY Fiyat: Bize Ulaşın
Teslimat tarihi için arayın

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

Teknik özellikler

Gövde malzemesi:	Glass
Hacim / Paket boyutu:	1 each
Kapasite:	25 mL
Ölçüm metodu:	To Contain
Sertifikalar:	Class A
Tür:	Tek Metrik, Karıştırma