



Bottle, storage, autoclavable, wide mouth, 125 mL pk/12

Ürün Numarası:

2939103

TRY Fiyat:

Bize Ulaşın

Mevcut

Polyethylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121°C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Teknik özellikler

| | |
|-----------------------|---------|
| Body Material: | Plastic |
| Hacim / Paket boyutu: | 12 /pk |
| Kapasite: | 125 mL |
| Malzeme: | PP |
| Package size UOM: | 12pk |