



## Flask, volumetric, class A, certified, 1000 mL

Ürün Numarası: 2636653

**ÜRETİMDEN KALDIRILAN ÜRÜN**  
Bu ürün artık mevcut değil.

Individually serialized glass flask supplied with Certificate of Graduation Accuracy and traceable to NIST Standards. Class A tolerances meet requirements of ASTM E438 and ASTM E288. Graduation ring on neck of flask. Supplied with ground glass stopper. (Stopper size 22).

### Teknik özellikler

|                  |               |
|------------------|---------------|
| Birim:           | 1 each        |
| Gövde malzemesi: | Class A Glass |
| Kapasite:        | 1000 mL       |
| Sertifikalar:    | Yes           |