



## Cylinder, graduated, mixing, certified, 50 mL

Ürün Numarası: 2636341  
TRY Fiyat: Bize Ulaşın  
Mevcut

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

### Teknik özellikler

Gövde malzemesi:	Glass
Hacim / Paket boyutu:	1 each
Kapasite:	50 mL
Ölçüm metodu:	To Contain
Sertifikalar:	Class A
Tür:	Tek Metrik, Karıştırma