



Disposable BOD bottles, 300 mL, 100/cs

Ürün Numarası: 2943100
TRY Fiyat: Bize Ulaşın
Mevcut

Disposable BOD Bottle

Case of 100 disposable BOD bottles **without stoppers**. Plastic bottle is EPA approved for BOD measurement.

Single bottles available under article number 2943101. Plastic stoppers (25 pack) available under article number 2943900. Glass stoppers (quantity of 1) available under article number 190901.

Fresh bottle, ready for use

No cleaning required

Unbreakable and recyclable

Oxygen will not diffuse in or out of sample

EPA compliant

Teknik özellikler

Body Material:	Plastic
Hacim / Paket boyutu:	100 100/pk
Kapasite:	300 mL
Malzeme:	PET
Metot numarası:	8043