



## Beaker, polypropylene, low form, 600 mL

Ürün Numarası: 108052  
TRY Fiyat: Bize Ulaşın  
Mevcut

Griffin low form. For general laboratory use. Non-toxic, noncontaminating, and biologically inert. Autoclavable. Combines no drip pouring, durability, and translucency. Approximate volume scales molded in for small sizes. Sizes 250 mL and larger feature silk-screened graduations and double spout for easier pouring. Nalgene.

### Teknik özellikler

Birim:	1 each
Çap:	110.7 mm with spout
Cihaz:	All except DR/4000, PC II
Kapasite:	600 mL
Malzeme:	Polypropylene, Low Form
Metot numarası:	8006