



pH paper, 6.0 - 8.0 pH range

Ürün Numarası:	37933
TRY Fiyat:	Bize Ulaşın
Mevcut	

Provides simple, fast, economical method of pH measurement. Sharp color change for each 0.3 pH unit. Convenient for sample pH adjustments. Recommended for sewage water, not finished water. Five 4.6-meter (15') rolls per package.

Teknik özellikler

Hacim / Paket boyutu:	5 4.6-meter (15') rolls per package
Ölçüm aralığı:	6.0 - 8.0 pH units