



Be Right™



pH paper, 1.0 - 11.0 pH range

Ürün Numarası:

39133

TRY Fiyat:

Bize Ulaşın

Mevcut

Provides simple, fast, economical method of pH measurement. Sharp color change for each 0.5 pH unit. Convenient for sample pH adjustments. Recommended for sewage water, not finished water. Five 4.6-meter (15') rolls per package.

Teknik özellikler

Hacim / Paket boyutu:

5 4.6-meter (15') rolls per package

Ölçüm aralığı:

1.0 - 11.0 pH units