



Be Right™



Gas washing bottle, 1200 mL

Ürün Numarası: 2662200
TRY Fiyat: Bize Ulaşın
Teslimat tarihi için arayın

For use in precipitating elements from gas. Used in the Hach Cold Vapor Mercury Apparatus Set (Prod. No. 2674400) to perform Cold Vapor Mercury Concentration Analysis–Hach Method 10065.

Teknik özellikler

Body Material:	Glass
Hacim / Paket boyutu:	1 each
Kapasite:	1200 mL
Malzeme:	Borosilicate Glass
Package size UOM:	1pk