



Fixed volume, Micropipette, 1000 µL

Ürün Numarası: 2946010
TRY Fiyat: Bize Ulaşın
Mevcut

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and color-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

Teknik özellikler

Doğruluk: 0,5 %
Hacim / Paket boyutu: 1 each
Kapasite: 1000 µL