



Test tube, rack, polyethylene, 16 mm (O.D.), 60 holes

Ürün Numarası: 2497901
TRY Fiyat: Bize Ulaşın
Mevcut

Sixty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5"). Color may vary.

Teknik özellikler

Çap:	16 mm
Hacim / Paket boyutu:	1 each
Test tüpü düzenlemesi:	5 x 12
Test tüpü yuvası:	60