



Rack, Test Tube, Polyethylene, 13 mm (O.D.), 90 Holes

Ürün Numarası:

2497900

TRY Fiyat:
Mevcut

Bize Ulaşın

Ninety-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Teknik özellikler

Çap:	13 mm
Hacim / Paket boyutu:	1 each
Test tüpü düzenlemesi:	6 x 15
Test tüpü yuvası:	90