



Test tube, rack, polyethylene, 25 mm (O.D.), 24 holes

Ürün Numarası: 2497903
TRY Fiyat: Bize Ulaşın
Mevcut

Twenty to ninety-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Teknik özellikler

Çap:	25 mm
Hacim / Paket boyutu:	1 each
Test tüpü düzenlemesi:	3 x 8
Test tüpü yuvası:	24