



Test kit, corrosion management, model CM-1

Ürün Numarası: 2495700

ÜRETİMDEN KALDIRILAN ÜRÜN
Bu ürün artık mevcut değil.

To reduce lead and copper levels in distribution, service and domestic water lines, suppliers must stabilize aggressive water by adjusting pH, temperature, and hardness. Analysts can verify the success of their stabilization by screening samples with this rugged, portable test kit.

Uses drop count titration and color disc methods

Maintain corrosion control of aggressive distribution system waters to meet Lead and Copper Rule Requirements

Establish baseline conditions to monitor treatment or disinfection strategy changes

Calculate Langlier and Agressive Indexes

Teknik özellikler

Platform:	Renk Diski
Test sayısı:	100
Uygulama:	Corrosion Management