



## Test kit, general boiler water, drop count titration, model AL-94 MG-L

Ürün Numarası:

222403

TRY Fiyat:  
2 hafta içinde teslimat

Bize Ulaşın

Tehlikeli



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Used for controlling boiler and cooling water treatment and checking corrosion, this test kit is an economical solution for the small boiler operator who needs to monitor alkalinity, chloride, hardness, iron, and pH.

**Uses simple, economical drop count titration and color disc methods**

**Kit contains all required reagents and apparatus to perform analysis in the field**

**Uses drop count titration and color disc methods**

---

### Teknik özellikler

Model:	AL-94 MG-L
Parametre:	Drop Count Titration
Platform:	Renk Diski
Uygulama:	Boiler Water