



## Test kit, professional boiler treatment/boiler feed, model BT-DT

Ürün Numarası: 2350600  
TRY Fiyat: Bize Ulaşın  
Teslimat tarihi için arayın

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. This kit provides testing capability for alkalinity, chloride, hardness, sulfite, phosphate, and pH, with a wide choice of ranges.

**Offers testing of treatment chemicals and feedwater parameters, with the accuracy demanded by water treatment professionals**

**The included Hach Digital Titrator simplifies titration procedures, saves time, and provides greater accuracy than drop count testing**

**Kit also includes color disc for simple, fast analysis**

---

### Teknik özellikler

Platform: Dijital Titratör  
Uygulama: Boiler Water