



Test kit, saltwater aquaculture, model FF-3

Ürün Numarası:

243003

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Maintaining good water quality is essential for profitable fish production. In many cases, the availability and cost of "crop" insurance is directly related to testing programs.

Designed to meet the demands of specialized salt-water applications

Each kit includes the apparatus and reagents necessary to measure 10 critical parameters using digital titrator, titration, and color disc methods

Two tests, ammonia-nitrogen and a convenient chloride titration test for salinity, have been designed exclusively for this kit

Teknik özellikler

Model:	FF-3
Platform:	Dijital Titratör
Uygulama:	Aquaculture