



Test kit, surface water

Ürün Numarası:

2559800

TRY Fiyat:

Bize Ulaşın

Teslimat tarihi için arayın

Tehlikeli



The ideal tool for source water monitoring and wastewater treatment plants, this kit uses drop count titration and color disc methods to provide economical, easy-to-use tests for ammonia, chlorine, pH, nitrate, dissolved oxygen, phosphorus, and temperature. The rugged, lightweight kit includes the convenient pH Pocket Pal Tester and all necessary apparatus and reagents for testing.

Ideal tool for both industrial and municipal water treatment plants

Uses drop count titration and color disc methods to provide economical, easy-to-use tests for ammonia, chlorine, pH, titrate, dissolved oxygen, phosphorus, and temperature

Rugged, lightweight kit includes the convenient pH Pocket Pal™ Tester and all necessary apparatus and reagents

Teknik özellikler

Platform:	Renk Diski
Test sayısı:	100
Uygulama:	Surface Water