



Chromium color disc test kit, model CH-8

Ürün Numarası:

183400

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Widely used as a corrosion inhibitor in cooling towers and engine cooling systems, chromium is also used in surface finishing, leather tanning, and other industries. The most widely used form, hexavalent chromium, is very toxic and strict limits are placed on discharges. Utilizing the color disc method to determine chromium (Range: 0.1-1.5 mg/L), this kit contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Metot:	Color disc/ Diphenylcarbazine
Model:	CH-8
Muhafaza Türü:	D
Ölçüm aralığı:	0 - 1,3 Cr ⁶⁺
Parameter:	Chromium Hexavalent
Parametre:	Chromium, low range - As hexavalent Cr
Platform:	Renk Diski
Smallest increment:	0.1
Test kiti modeli:	CH-8
Test sayısı:	100