



Detergents test kit, model DE-2

Ürün Numarası:

143203

TRY Fiyat:

Bize Ulaşın

Teslimat tarihi için arayın

Tehlikeli



The presence of detergents in surface or ground waters indicates industrial or domestic pollution. In natural waters anionic LAS and ABS detergents cause excessive foaming, produce undesirable taste and odors, promote growth of disease-producing organisms, and tend to kill aquatic life. This kit's simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Teknik özellikler

Metot:	Color disc/
	Toluidine blue-O Chloroform
Model:	DE-2
Muhafaza Türü:	Q
Ölçüm aralığı:	0 - 1,2 mg/L LAS
Parametre:	Detergents - As LAS and/or ← ABS detergents
Platform:	Renk Diski
Test kiti modeli:	DE-2
Test sayısı:	32