



Test kit, phosphonate, model PN-10, w/out uv lamp & power supply

Ürün Numarası:

2113301

TRY Fiyat:
3 gün içinde teslimat

Bize Ulaşın

Tehlikeli



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Contains 100 tests.

Simple, fast, color disc method

UV lamp digestion converts phosphonates into reactive phosphate

Control phosphonate concentration to save chemical costs

Hach's reliable PhosVer 3 reagent provides results in under 15 minutes

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Güç kaynağı:	No UV Lamp Nor Power Supply
Metot:	Color disc
Model:	PN-10
Muhafaza Türü:	F
Ölçüm aralığı:	0 - 100 mg/L PO ₄
Parameter:	Phosphonate
Parametre:	Phosphonate - As PO ₄ ³⁻
Platform:	Renk Diski
Smallest increment:	1
Standart aralık:	0 - 100 mg/L PO ₄
Test kiti modeli:	PN-10
	(As above but without UV lamp and power supply.)
Test sayısı:	0 - 100