



Iron (Ferrous) Color Disc Test Kit, Model IR-18C, 0-7 mg/L

Ürün Numarası:

2667200

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



The most accurate visual test kit method, Hach color discs provide fast, accurate comparisons. Range: 0-7 mg/L. Contains 100 tests.

Simple, fast, color disc method

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Teknik özellikler

Metot:	Color disc/ 1,10 Phenanthroline
Model:	IR-18C
Muhafaza Türü:	B
Ölçüm aralığı:	High Single Range
Parameter:	Iron Ferrous
Parametre:	Iron, medium range - As Fe ²⁺ (ferrous iron)
Platform:	Renk Diski
Smallest increment:	0.2
Test kiti modeli:	IR-18C
Test sayısı:	100
Tür:	Tek Metrik