



Manganese test kit, model MN-5

Ürün Numarası:

146700

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. An economical solution, this kit uses the simple, accurate color disc method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Teknik özellikler

Metot:	Color disc/Cold periodate
Model:	MN-5
Muhafaza Türü:	D
Ölçüm aralığı:	0 - 3,0 mg/L Mn
Parameter:	Manganese
Parametre:	Manganese - As Mn
Platform:	Renk Diski
Smallest increment:	0.1
Standart aralık:	0 - 3.0 mg/L Mn
Test kiti modeli:	MN-5
Test sayısı:	100