



Test kit tannin/lignin

Ürün Numarası:

193701

TRY Fiyat:

Bize Ulaşın

3 gün içinde teslimat

Tehlikeli



A product of decomposed plant material, tannin occurs in natural waters and is used for boiler and cooling tower water treatment, dyeing, ink manufacture, tanning, and paper sizing. Lignin, also a plant product, is found in paper pulp manufacturing waste. Tannin and lignin are analyzed together and reported as mg/L tannic acid. This kit utilizes the simple, low-cost color disc method - the most accurate visual test kit method.

Tannin and lignin are analyzed together and reported as mg/L tannic acid.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Teknik özellikler

Metot:	Color disc/Tannin/Lignin
Model:	Tannin/Lignin
Muhafaza Türü:	R
Ölçüm aralığı:	0 - 15 mg/L as Tannic Acid
Ölçüm aralığı (2):	0 - 150 mg/L
Parameter:	Tannin/Lignin
Parametre:	Tannin/Lignin - As tannic acid
Platform:	Renk Diski
Smallest increment:	0.5
	5
Standart aralık:	0 - 15 mg/L as Tannic Acid; 0 - 150 mg/L
Test kiti modeli:	TA-3
Test sayısı:	100