



Test kit, oxygen scavenger, color disc, model OS-1, 100 tests

Ürün Numarası:

2349500

Tehlikeli

TRY Fiyat:

Bize Ulaşın

2 hafta içinde teslimat



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 100 tests.

Simple and fast color disc method

Measures common oxygen scavengers: DEHA, hydroquinone, ISA, and MEKO

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Metot:	Color disc/Iron reduction
Model:	OS-1
Muhafaza Türü:	R
Ölçüm aralığı:	0 - 52 µg/L DEHA (and other Oxygen Scavengers)
Ölçüm aralığı (2):	0 - 1360 µg/L DEHA (and other Oxygen Scavengers)
Ölçüm aralığı (3):	0 - 300 µg/L DEHA (and other Oxygen Scavengers)
Parametre:	Oxygen Scavenger - As DEHA, hydroquinone, erythorbic acid, and methylethylketoxime
Platform:	Renk Diski
Test kiti modeli:	OS-1
Test sayısı:	50