



Test kit, phenols, model PL-1

Ürün Numarası:

2483600

TRY Fiyat:

Bize Ulaşın

2 hafta içinde teslimat

Tehlikeli



Utilizing the simple color disc method with a range of 0 - 4 mg/L, this kit contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Metot:	Color disc/ 4-aminoantipyrine
Model:	PL-1
Muhafaza Türü:	D
Ölçüm aralığı:	0 - 4,0 mg/L
Parameter:	Phenols
Parametre:	Phenols
Platform:	Renk Diski
Smallest increment:	0.2 mg/L
Standart aralık:	0 - 4.0 mg/L
Test kiti modeli:	PL-1
Test sayısı:	100