



## Phosphorus, orthophosphate (reactive) test kit, model PO-14

Ürün Numarası:

147500

Tehlikeli

TRY Fiyat:  
Mevcut

Bize Ulaşın



Simple, low-cost color disc method is the most accurate visual test kit method

### Teknik özellikler

|                     |                                                        |
|---------------------|--------------------------------------------------------|
| Metot:              | Color disc/Stannous                                    |
| Model:              | PO-14                                                  |
| Muhafaza Türü:      | D                                                      |
| Ölçüm aralığı:      | Single Range                                           |
| Parameter:          | Phosphorus                                             |
| Parametre:          | Phosphorus, Orthophosphate (reactive) - As $PO_4^{3-}$ |
| Platform:           | Renk Diski                                             |
| Smallest increment: | 0.1                                                    |
|                     | 1                                                      |
| Standart aralık:    | 0 - 4.4 mg/L $PO_4$                                    |
| Test kiti modeli:   | PO-14                                                  |
| Test sayısı:        | 100 C**                                                |