



Test kit color disc chlorine dioxide, model CLD-2

Ürün Numarası:

2842800

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Monitor ClO_2 levels on-site with the NEW Chlorine Dioxide Test Kit. The Chlorine Dioxide Test Kit uses the reliable and EPA-accepted DPD method to monitor ClO_2 concentration, while eliminating chlorine interference through the addition of glycine. With 100 tests per kit and available range of 0-2 mg/L and resolution of 0.1 mg/L, operators can accurately monitor around the 0.8 mg/L regulated limit.

Simple, cost-effective method for maintaining optimum chlorine dioxide levels in your plant

USEPA-accepted for waste water and potable water

Approved by some states for potable water under specific conditions (check with your regulatory agency)

Teknik özellikler

Metot:	Color Disc/DPD
Model:	CLD-2
Muhafaza Türü:	D
Ölçüm aralığı:	0 - 2,0 ClO_2
Parameter:	Chlorine Dioxide
Parametre:	Chlorine Dioxide
Platform:	Renk Diski
Smallest increment:	0.1
Standart aralık:	0 - 2.0 ppm ClO_2
Test kiti modeli:	CLD-2 ²
Test sayısı:	100