



Nitrate test kit, model NI-11

Ürün Numarası:

146803

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



One of Hach's most popular test kits, it uses the simple, low-cost color disc method.

One of Hach's most popular test kits

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Teknik özellikler

Metot:	Color disc/ Cadmium reduction
Model:	NI-11
Muhafaza Türü:	B
Ölçüm aralığı:	Single Range
Parameter:	Nitrogen, Nitrate
Parametre:	Nitrogen, Nitrate - As NO ₃ ⁻ -N
Platform:	Renk Diski
Smallest increment:	1.0
Standart aralık:	0 - 40 mg/L NO ₃ ⁻ -N
Test kiti modeli:	NI-11**
Test sayısı:	100