



Test kit, nitrogen , model NI-8

Ürün Numarası:

224100

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Designed for testing ammonia in freshwater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

One of Hach's most popular test kits

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Teknik özellikler

| | |
|---------------------|-----------------------------------------------|
| Metot: | Color disc/Nessler reagent |
| Muhafaza Türü: | B |
| Ölçüm aralığı: | 0 - 3,0 mg/L NH ₃ -N |
| Parametre: | Nitrogen, Ammonia, mid range - For freshwater |
| Platform: | Renk Diski |
| Smallest increment: | 0.1 |
| Test kiti modeli: | NI-8 |
| Test sayısı: | 100 |