



Test kit, nitrogen, ammonia, model NI-SA

Ürün Numarası:

2428700

TRY Fiyat:

Bize Ulaşın

3 gün içinde teslimat

Tehlikeli



Designed for testing ammonia in freshwater or seawater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

Designed for testing ammonia in freshwater or seawater

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Teknik özellikler

Metot:	Color disc/Salicylate
Muhafaza Türü:	D
Ölçüm aralığı:	0 - 2,4 mg/L NH ₃ -N
Parametre:	Nitrogen, Ammonia, mid range - For freshwater or seawater**
Platform:	Renk Diski
Smallest increment:	0.1
Test kiti modeli:	NI - SA
Test sayısı:	100