



Test kit, phosphonate, model PN-10, w/ 230 vac uv lamp & power supply

Ürün Numarası:

2113302

TRY Fiyat:

Bize Ulaşın

3 gün içinde teslimat

Tehlikeli



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Contains 100 tests.

Simple, fast, color disc method

UV lamp digestion converts phosphonates into reactive phosphate

Control phosphonate concentration to save chemical costs

Hach's reliable PhosVer 3 reagent provides results in under 15 minutes

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Güç kaynağı:	230 Vac UV Lamp and Power Supply
Metot:	Color disc/UV photolysis
	(230 Vac, 50 Hz)
Model:	PN-10
Muhafaza Türü:	F
Ölçüm aralığı:	0 - 100 mg/L PO_4
Parametre:	Phosphonate - As PO_4^{3-}
Platform:	Renk Diski
Smallest increment:	1
Test kiti modeli:	PN-10
Test sayısı:	0 - 100