



## Test kit, molybdate, model MO-LR

Ürün Numarası:

2359300

TRY Fiyat:  
Mevcut

Bize Ulaşın

Tehlikeli



Molybdates are gaining over chromates for use as corrosion inhibitors due to lower toxicity. The high cost of these treatments warrant careful control of their concentration. Kit reagents react with Molybdate to form a yellow color, which is measured with a color disc.

**Kit reagents react with Molybdate to form a yellow color, which is measured with an easy-to-use color disc method**

**Effectively monitors corrosion inhibitors**

### Teknik özellikler

Metot:	Color disc/ Ternary complex
Model:	MO-LR
Muhafaza Türü:	D
Ölçüm aralığı:	Single Range
Parameter:	Molybdate
Parametre:	Molybdate - As Molybdenum
Platform:	Renk Diski
Smallest increment:	0.25
Standart aralık:	0 - 3.0 mg/L Mo
Test kiti modeli:	MO-LR
Test sayısı:	100