



Test Kit, Accuvac Color Disc, Ozone, 0.05 - 1.50 mg/L

Ürün Numarası:

2518050

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



This ozone test kit uses the fast, accurate visual method for determining ozone concentrations. Use to control disinfectant levels in water, wastewater, and industrial water. Kit contains 25 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

AccuVac vials mean no mess, no measuring

Indigo trisulfonate method masks interferences from chlorine and other oxidizing agents and makes ozone analysis easier

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Teknik özellikler

Metot:	Color disc/ Indigo trisulfonate
Model:	NA
Muhafaza Türü:	G
Ölçüm aralığı:	0 - 1,50 mg/L O ₃
Parameter:	Ozone
Parametre:	Ozone - As O ₃
Platform:	Accuvac
Smallest increment:	0.05
Standart aralık:	0 - 1.50 mg/L O ₃
Test kiti modeli:	AccuVac® Kit
Test sayısı:	25