



Test kit, copper (free) test color disc, model CU-5

Ürün Numarası:

1421300

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Copper may exist in natural waters and effluents as a soluble salt or as suspended solids. A small amount is essential for plants and animals, and concentrations exceeding 0.1 mg/L are useful for controlling algae and plankton growth. Quantities ranging from 0.02 to 0.1 mg/L are toxic for some fish, so its use for treating fish disease requires careful monitoring. This kit is USEPA-accepted for waste water and potable water and approved by some states for potable water under specific conditions (Check with your regulatory agency). Range: 0.1-5 mg/L; contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Metot:	Color-disc/Bicinchoninate
Model:	CU-5
Muhafaza Türü:	B
Ölçüm aralığı:	0 - 4 mg/L Cu
Parameter:	Copper (free)
Parametre:	Copper - As free Cu
Platform:	Renk Diski
Smallest increment:	0.1
Standart aralık:	0 - 4 mg/L Cu
Test kiti modeli:	CU-5
Test sayısı:	100 F