



## Test kit, copper (free & total) color disc, model CU-6

Ürün Numarası:

2194100

TRY Fiyat:  
Mevcut

Bize Ulaşın

Tehlikeli



Copper may exist in natural waters and effluents as a soluble salt or as suspended solids. A small amount is essential for plants and animals, and concentrations exceeding 0.1 mg/L are useful for controlling algae and plankton growth. Quantities ranging from 0.02 to 0.1 mg/L are toxic for some fish, so its use for treating fish disease requires careful monitoring. Range: 0.1-4 mg/L; Contains 50 Free and 50 Total tests.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Teknik özellikler

Metot:	Color-disc/Bicinchoninate, Hydrosulfite reduction
Model:	CU-6
Muhafaza Türü:	B
Ölçüm aralığı:	0 - 4 mg/L Cu
Parameter:	Copper (free & total)
Parametre:	Copper - As free & total Cu <sup>#</sup>
Platform:	Renk Diski
Smallest increment:	0.1
Test kiti modeli:	CU-6
Test sayısı:	50 (Free), 50 (Total)