



Test kit, formaldehyde, model FM-1

Ürün Numarası:

2183100

TRY Fiyat:

Bize Ulaşın

Teslimat tarihi için arayın

Tehlikeli



Formaldehyde solutions are used for maintenance of reverse osmosis systems and as disinfectants. The Model FM-1 measures higher concentrations by drop count titration. Contains 100 tests.

Utilises drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Metot:	Drop count titration/ Thymolphthalein
Muhafaza Türü:	D
Ölçüm aralığı:	0,05 - 1% & 0,5 - 10%
Parametre:	Formaldehyde - As CH_2O
Platform:	Drop Count Titration
Smallest increment:	0.05% 0.5%
Test kiti modeli:	FM-1
Test sayısı:	100