



Test kit, hydrogenperoxide, model HYP-1

Ürün Numarası:

2291700

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Model HYP-1 Hydrogen Peroxide Test Kit (as H_2O_2) utilizes drop count titration and contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Metot:	Drop count titration/ Thiosulfate
Muhafaza Türü:	D
Ölçüm aralığı:	0,2 - 2,0 mg/L & 1 - 10 mg/L
Parametre:	Hydrogen Peroxide - As H_2O_2
Platform:	Drop Count Titration
Smallest increment:	0.2
	1
Test kiti modeli:	HYP-1
Test sayısı:	100