



Test kit, chlorine (total), model CN-65

Ürün Numarası:

225401

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

En küçük artış adımları:

1

Metot:

Drop count titration/

Thiosulfate

Muhafaza Türü:

D

Ölçüm aralığı:

0,2 - 4,0 mg/L & 1 - 20 mg/L

Parametre:

Chlorine, mid and high range - As Total Cl₂

Platform:

Drop Count Titration

Smallest increment:

0.2

1

Test kiti modeli:

CN-65

Test sayısı:

100