



Test kit, hypochlorite, model CN-HR

Ürün Numarası:

2687200

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Hypochlorite Kit, Model CN-HR - Drop Count

Utilizing the simple, economical drop count titration method, this kit is used to determine high range hypochlorite solution concentrations (bleach solutions). Range: 5-15%, contains 100 tests.

Simple, economical titration method

Kit contains all required reagents, solutions, and glassware

Teknik özellikler

Metot:	Drop count titration/ Thiosulfate
Muhafaza Türü:	D
Ölçüm aralığı:	50 - 150 g/L Cl_2 (5-15%)
Parametre:	Hypochlorite, high range - (Bleach) As Cl_2
Platform:	Drop Count Titration
Smallest increment:	0.5%
Test kiti modeli:	CN-HR
Test sayısı:	100