



## Test kit, chromium, model CH-14

Ürün Numarası:

222702

TRY Fiyat:

Bize Ulaşın

2 hafta içinde teslimat

Tehlikeli



Widely used as a corrosion inhibitor in cooling towers and engine cooling systems, chromium is also used in surface finishing, leather tanning, and other industries. The most widely used form, hexavalent chromium, is very toxic and strict limits are placed on discharges. Utilizing the drop count titration method to determine chromium (high range), this kit contains 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Teknik özellikler

Metot:	Drop count titration/ Thiosulfate
Muhafaza Türü:	D
Ölçüm aralığı:	5 - 100 mg/L & 50 - 1000 mg/L Cr <sup>6+</sup>
Parametre:	Chromium, high range - As Hexavalent Cr
Platform:	Drop Count Titration
Smallest increment:	5 50
Test kiti modeli:	CH-14
Test sayısı:	100