



Test kit, hardness (total), model HA-DT

Ürün Numarası:

2063600

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Utilizing the Hach Digital Titrator, the Model HA-DT kit is used to determine hardness concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-4,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Teknik özellikler

En küçük artış adımları:	0.1-10
Metot:	Digital Titrator/EDTA
Muhafaza Türü:	F
Ölçüm aralığı:	10 - 4,000 CaCO ₃
Parametre:	Hardness, Total - As CaCO ₃
Platform:	Dijital Titratör
Test kiti modeli:	HA-DT
Test sayısı:	100