



## Test kit, chloride, model CD-DT

Ürün Numarası:  
TRY Fiyat:  
Teslimat tarihi için arayın

2063500  
Bize Ulaşın

Tehlikeli



The single-parameter Model CD-DT uses the Hach Digital Titrator with Mercuric Nitrate Titrant Cartridges to determine chloride concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-8,000 mg/L. Contains approximately 100 tests.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

### Teknik özellikler

Metot:	Dijital Titrator/ Mercuric nitrate
Muhafaza Türü:	F
Ölçüm aralığı:	10 - 8,000 mg/L Cl <sup>-</sup>
Parametre:	Chloride - As Cl <sup>-</sup> , 10-8,000 mg/L Cl <sup>-</sup>
Platform:	Dijital Titratör
Smallest increment:	0.1-20
Test kiti modeli:	CD-DT
Test sayısı:	100